



- Expressly designed for drying and curing of soft gels and capsules
- Dropped sides and/or ends improve air circulation to speed drying
- Easy stacking without warping or distortion
- Smooth, non-porous and easy to clean
- Autoclavable material options
- Temperature range continuous -60° to 250° F (-51° to 121° C)
- Heavy-duty dollies available

## MFG Gel Cap Drying Trays

Model No.	Description	Top Outside Dimensions (L x W x H)					
		inches	mm	inches	mm	inches	mm
600008	Tray with drop ends	29 7/8	759	23 7/8	606	1 1/2	38
600108	Tray with drop sides	29 7/8	759	23 7/8	606	1 1/2	51
600208	Tray with drop sides	23 7/8	606	14 15/16	379	1 3/8	35
600308	Tray with drop ends	29 7/8	759	23 7/8	606	2	51
600408	Tray with drop ends	29 7/8	759	23 7/8	606	1	25
642008	Tray with drop ends and drop sides	29 7/8	759	23 7/8	606	1 3/8	35
643008	Tray with drop sides	32 7/8	833	18 7/8	478	1 5/8	39

**Standard Colors: White, Gray, Beige, Light Blue, Light Green**

### Standard Colors



Gel Cap Drying Trays  
Model No. 600 Stack



Gel Cap Drying Tray  
Model No. 642008

Please contact your preferred distributor to order.



MOLDED FIBER GLASS TRAY  
6175 Route 6 · Linesville, PA 16424 USA

PH 814 683.4500 or 800 458.6050  
FX 814 683.4504    www.mfgtray.com

MFGTPHgcdt09/23

